RADIANT WARMER

Instructions to turn on machine and turn off radiant heater.

HOW TO POWER ON AND PREVENT ALARMING

Step 1.) Press button on rear of display to turn on warmer.

Step 2.) Press button on front of display to turn on light.

Step 3.) On screen, under manual press and hold DOWN ARROW to remove all green bars.

Step 4.) Press OK button to continue.

Step 5.) On screen, under mattress °C press and hold DOWN ARROW lowering number to 36.0.

Step 6.) Press OK button to continue.

Step 7.) Again, under mattress °C press and hold

DOWN ARROW to lower number to 30.0.

Step 8.) Press OK button.

Click to HERE to view item details on Alpha Props website

