

HRD EKG Monitor

How to Boot into DEMO MODE

A video walk-through of these instructions can be viewed at <https://youtu.be/vU7XB3NupdY>

Step 1.) Switch on power.

Step 2.) Press "Main" button

A menu will open.

Step 3.) Dial knob to move through menu to "Demo" and press knob to select.

Step 4.) Press knob a second time to highlight password field.

Step 5.) Dial knob to change number to "5188" and press to select.

You are now in DEMO MODE.

[Click to HERE](#) to view item details on
Alpha Props website

